

# Restriction Enzyme BsrF I



Cat.# FG-BsrFI

Size 1,000 units Conc. 10 units/µl

Store at -20℃

Supplied with: 10X FastGene® Buffer IV (FG-REB4)

10X FastGene® FastCut Buffer (FG-REBHF) 6X DNA Loading Buffer

Sterile water

### Recognition site

For Research Use Only. Not for use in diagnostic procedures.

**ISO**9001

Source: Bacillus stearothermophilus CPW16

#### Reaction conditions

1X FastGene® Buffer IV 37°C 1X FastGene® FastCut Buffer, 37°C

## FastGene® FastCut Buffer

FastGene® restriction enzyme can cut substrate DNA in 5-15 with FastGene® FastCut Buffer.

#### 1X FastGene® Buffer IV

20 mM Tris-acetate (pH 7.9 at 25°C) 50 mM potassium acetate 10 mM magnesium acetate 100 μg/ml BSA

#### Unit definition

One unit is defined as the amount of enzyme required for complete digestion of 1  $\mu g$  pBR322 at 37°C for 1 hr in 50  $\mu l$  reaction mixtures.

## Quality control

- Unit definition assay
- Overdigestion assay
- Endonuclease assay
- Extreme pure assay

# Dilution buffer:

FastGene® Diluent C

## **Heat Inactivation**

No.

#### Methylation sensitivity

dam methylation: Not sensitive dcm methylation: Not sensitive CpG methylation: sensitive

## Prolonged incubation

A minimum amount of enzyme required to digest 1  $\mu$ g substrate DNA for 16 hr; 0.13 U.

# Relative activity in FastGene® Buffers

 FastGene® Buffer I:
 75%

 FastGene® Buffer II:
 100%

 FastGene® Buffer III:
 100%

 FastGene® Buffer IV:
 100%

 FastGene® FastCut Buffer:
 100%

#### Note

It is an isoschizomer of Cfr10 I. Cleavage of mammalian genomic DNA is blocked by CpG methylation. Reaction condition of low salt, excess enzyme, excess glycerol (>5%) or high pH (>8.0) may result in star activity.

# Standard reaction condition

- Normal protocol

Component	Final Conc.	Volume
Substrate DNA	1 μg	ΧμΙ
10X FastGene® Buffer IV	1 X	5 μΙ
BsrF I	10 unit	1 μΙ
Sterile water		up to 50 μl

→ Incubate at 37°C for 1 hr

- Fast protocol

Component	Final Conc.	Volume
Substrate DNA	1 μg	ΧμΙ
10X FastGene® FastCut Buffer	1 X	5 μΙ
BsrF I	10 unit	1 μΙ
Sterile water		up to 50 μl

 $\rightarrow$  Incubate at 37°C for 15 min