

# Restriction Enzyme Eco47 I



Cat.# FG-Eco47I Size 800 units Conc. 10 units/µl

Store at -20℃

Supplied with: 10X FastGene® Buffer III (FG-REB3)

10X FastGene® FastCut Buffer (FG-REBHF)

6X DNA Loading Buffer

Sterile water

## Recognition site

For Research Use Only. Not for use in diagnostic procedures.

**ISO**9001

#### Dilution buffer:

FastGene® Diluent A

#### Heat Inactivation

Eco47 I can be inactivated at 65°C for 20 min.

## Methylation sensitivity

dam methylation: Not sensitive dcm methylation: sensitive CpG methylation: sensitive

# Relative activity in FastGene® Buffers

 FastGene® Buffer I:
 100%

 FastGene® Buffer II:
 100%

 FastGene® Buffer III:
 100%

 FastGene® Buffer IV:
 100%

 FastGene® FastCut Buffer:
 100%

Source: E.coli RFL47

## Reaction conditions

1X FastGene® Buffer III 37°C 1X FastGene® FastCut Buffer, 37°C

## FastGene® FastCut Buffer

FastGene® restriction enzyme can cut substrate DNA in 5-15 with FastGene® FastCut Buffer.

# 1X FastGene® Buffer III

50 mM Tris-HCl (pH 7.9 at 25°C) 100 mM NaCl 10 mM MgCl $_2$  100  $\mu$ g/ml BSA

#### Unit definition

One unit is defined as the amount of enzyme required for complete digestion of 1  $\mu$ g bacteriophage  $\lambda$  at 37°C for 1 hr in 50  $\mu$ l reaction mixtures.

# Quality control

- Unit definition assay
- Overdigestion assay
- Endonuclease assay
- Extreme pure assay

#### Standard reaction condition

- Normal protocol

Component	Final Conc.	Volume
Substrate DNA	1 μg	Xμl
10X FastGene® Buffer III	1 X	5 µl
Eco47 I	10 unit	1 µl
Sterile water		up to 50 μl

→ Incubate at 37°C for 1 hr

- Fast protocol

Component	Final Conc.	Volume
Substrate DNA	1 μg	Χ μΙ
10X FastGene® FastCut Buffer	1 X	5 μΙ
Eco47 I	10 unit	1 μΙ
Sterile water		up to 50 μl

→ Incubate at 37°C for 15 min

 $\ensuremath{\mathbb{X}}$  We recommend 5-10 units of enzyme per  $\mu g$  DNA and 10-20 units for genomic DNA in a 1 h digest.